

DNA Extraction

Fruit: Kiwi fruit
Strawberries.
Bananas

Equipment: Mortar and pestle.
Knife. scissors
250 ml. beaker.
Glass rod.
Salt
Large test tubes x2 / stopper
Filter funnel / filter paper.
50 ml. measuring cylinder
Methylated spirits (35 mls.)

Method:

1. Peel kiwi fruit. Cut into 4 parts with a knife.

(1 part/group)

Cut into small pieces with scissors.

Pulp using mortar and pestle.

Place kiwi fruit pulp into a large test tube. Set aside.

2. Add 200mls. of hot tap water (50°C) to a beaker and add 1 teaspoon of salt to it. Mix with a glass rod.

3. Add enough of the salty water to just cover the kiwi pulp in the large test tube.

Using a stopper in the test tube, shake for 20 secs.

4. Filter the shaken solution into another large test tube.

5. Slowly add about 35 mls. of cooled methylated spirits down the side of the test tube so as not to mix the meths./ extract layers.