


MICROSCOPES

The following diagram of the common microscope and its labelled parts need to be remembered.



Important Rules For Using Microscopes

- Make sure the mirror reflects light up through the stage.
- Focus on **low power first**.
- **Look from the side** to adjust the objective lens **just above the slide**.
- Use the **coarse focus** to wind the eye piece upwards under **low power**.
- **Look from the side when switching to higher power** lenses to avoid hitting the slide.
- Use the **fine focus** when looking under **high power**.

Activity

- Do Activity 3.1 p. 42 Science Quest 8.

Magnification

The two lenses that determine the **magnification** of your microscope are the eyepiece lens and the objective lens. Each lens has a number on it that signifies its magnification. Multiplying the eyepiece number by the objective lens number will give you the magnification of the microscope. For example:

- eyepiece: $\times 10$
- objective: $\times 40$
- magnification = $\times 400$.



Field of view 4 mm
(4000 μm)
Magnification $\times 40$



Field of view 1.6 mm
(1600 μm)
Magnification $\times 100$



Field of view 0.4 mm
(400 μm)
Magnification $\times 400$

As the field of view gets smaller, the magnification gets larger.