


RANGE OF RADIATION IN AIR

Aim: To determine how far alpha and beta particles penetrate through air.

Hypothesis:

Method:

1. Set up the following apparatus, without the source in place.




2. Place the Americium-241 source (alpha source) under the Geiger counter and measure the count rate for 60 seconds at the following distances: 1.0 cm, 2.0 cm, 3.0 cm, 4.0 cm and 5.0 cm. Record your observations in a table.
3. Repeat step 2 for the Strontium-90 source (beta source). Record your observations in a table.
4. Graph the results for both sources - use one graph if possible.

Results:

Alpha Source	
Distance (cm)	Count Rate (per minute)
1.0	
2.0	
3.0	
4.0	
5.0	

Beta Source	
Distance (cm)	Count Rate (per minute)
1.0	
2.0	
3.0	
4.0	
5.0	



Discussion:

Conclusion: